
Cube Mate™ Humidifier

Item #: 40-682-000

Instructions for Use: Please read these instructions carefully and keep them for future reference.

Includes: 1 Cube Mate™ Ultrasonic Cool Mist Humidifier

KEY TO SYMBOLS

- | | |
|--------------------------------------|-----------------------------|
| 1. Mist Output Opening | 8. Ultrasonic Transducer |
| 2. Aromatherapy Tray | 9. Air Inlet |
| 3. Water Tank | 10. Power Adapter Input |
| 4. Tank Cover | 11. Power Adapter |
| 5. Finger Holds | 12. Minimum Fill Level Line |
| 6. Power Knob / Adjustable Mist Knob | 13. Maximum Fill Level Line |
| 7. Power / Refill Indicator Light | 14. Corner Drainage Columns |

Important Safety Instructions:

- After removing all packaging materials, ensure that the humidifier has not been damaged during transportation.
- Do not leave any packaging materials within the reach of children to avoid injury.
- Do not move, tilt or tip the humidifier during use. Doing so may cause water leakage.
- Before connecting the equipment, ensure that the rated voltage (AC 100-240V) corresponds to the main voltage supply of the electrical outlet.
- Do not use additional power strip adapters, and/or extension cords to avoid the risk of fire, shock or injury.
- Only place the humidifier on flat and stable surfaces.
- Only use clean, cool tap, filtered or distilled water to fill the water tank. Do not use hot water. Do not add salt, fragrance/essential oils, medication or descaling tablets to the water or water tank.
- Do not operate on wood or other water sensitive surfaces, as condensation may occur.
- Do not leave the humidifier unattended in a closed room during use, as the air could become too saturated and leave condensation on furniture or walls.
- Always disconnect the equipment from the power outlet when it is left unattended or not in use.
- Before performing any cleaning or maintenance, disconnect the equipment from the power adapter and outlet.
- This product is not a toy. Keep out of reach of children.
- This equipment is a household humidifier for indoor use only. Do not operate outdoors.
- Do not use accessories, spare parts, or any component not supplied or approved by the manufacturer.
- Never touch the equipment with wet or damp hands.
- If the equipment is damaged and/or not working properly, turn the unit off by disconnecting the unit from the power adapter and the plug from the power outlet.
- Never submerge the equipment in water or any other liquid.
- Do not touch the equipment if it accidentally falls into the water.

- Do not direct the vapor flow towards furniture or electrical appliances.
- Do not touch the ultrasonic transducer with your hands or tools while the unit is plugged in.
- Do not operate the humidifier without water.
- Do not attempt to empty or fill the tank while the unit is plugged in.
- Clean and dry the humidifier at the end of every operating cycle.
- Dispose of any worn accessories and/or the equipment at an appropriate collection point, as required by state regulations.

How to prepare and use the equipment:

Check that the equipment is completely clean, the power knob is in the OFF position and that the power adapter is not connected to the unit or power outlet.

Gently remove the tank cover from the base using the finger holds, and fill the water tank to the peak level in between the minimum and maximum fill lines. Do NOT fill water past the maximum fill level line.

If fragrance is desired, remove the aromatherapy tray and place one or two drops of fragrance or essential oil directly onto the pad inside the tray. Replace the aromatherapy tray back into the top of the water tank.

Note: Do NOT add fragrance or essential oil directly to the water or water tank, as this may cause damage to the internal electrical components.

Replace the tank cover, making sure it is properly secured in place.

Connect the equipment to the adapter and connect to a power outlet.

Use the mist regulation knob to turn the unit ON by moving the knob clockwise, and the power indicator light will illuminate blue.

Vapor will start to exit.

Use the mist regulation knob to control the amount of vapor flow desired.

To turn the unit OFF, turn the knob counterclockwise.

When water has reached the minimum fill level in the tank, the humidifier stops automatically and the power refill indicator light will illuminate red.

Cleaning and Maintenance:

- Before cleaning the equipment, always turn the humidifier OFF, and disconnect the unit from the adapter and plug from the power outlet.
- Do not place the humidifier directly under running water.
- Do not wet the adapter or air inlet.
- Do not scrape the ultrasonic transducer with metallic materials, tools or objects.
- Do not leave water in the water tank when the equipment is not in use.

Daily Cleaning:

Remove the tank cover from the base and slowly pour out all remaining water.

Clean the cover, tank and transducer of the humidifier with a soft, clean cloth and dry thoroughly.

Dry all parts thoroughly before storage.

Weekly Cleaning:

Remove the tank cover from the base, and quickly pour out all remaining water.

Note: The unit is equipped with four drainage columns located on each inside corner of the unit. When pouring out the water, it may flow into these columns, thus causing water to drain from the bottom of the unit. It is recommended to always drain the water over a sink to avoid potential spills.

Pour a water solution prepared with 50% white vinegar into the water tank. Do not fill the tank all the way to the brim.

Replace the tank cover and leave soak for at least 4 hours (overnight is ideal). Do NOT turn the humidifier ON and run the unit with the water and vinegar cleaning solution.

Remove the tank cover and pour out all of the water/white vinegar solution, and rinse the tank with cool, clean water.

Wipe the cover and tank of the humidifier with a soft, clean cloth and dry thoroughly.

Clean the outer surfaces of the humidifier with a soft, damp cloth.

Do not use abrasive products, alcohol, solvents or similar products which may cause permanent damage to the outer surface.

To clean the aromatherapy tray, remove the pad and gently rinse thoroughly with water. Allow the pad to air dry before replacing into the tray.

Troubleshooting:

Problem	Solution
If the humidifier does not work	<ul style="list-style-type: none">• Check that the adapter and plug are connected correctly.• Check that the switch is in the ON position.• Check that there is sufficient water in the water tank (water fill level must be in between the minimum and maximum fill line).
If the power knob is in the ON position but there is no vapor flow	<ul style="list-style-type: none">• Check whether there are hard water deposits on the ultrasonic transducer.• Check whether there is sufficient water in the water tank.
Strange odor	<ul style="list-style-type: none">• Clean the water tank and reservoir with water and vinegar solution, then rinse with clean water.• Aerate the water tank in a cool, dry place for 12 hours.

Technical Specifications:

Run Time: Min. up to 5 hours/Max. up to 2 ½ hours

Tank Capacity: 8.45 fl. oz. between minimum and maximum fill lines

Ideal Room Size: Small

Average Humidity Output*: 55 ml/hr

Power Consumption: AC 100-240V, 50-60Hz, 12W

Sound Level: <35 dBA

*Humidity output may vary based on indoor room temperature and humidity, and the mist output setting.

This product has been tested and complies with the requirements for Federal Communications Commission, Part 18 for Radio/TV/Communication interference. Although tested, it may affect these devices. If the humidifier is found to interfere, separate the device and/or the humidifier. Conduct only the user maintenance found in this manual. Other maintenance and servicing can cause harmful interference and can void the required FCC compliance.

Intertek

Conforms to UL STD. 998.

Limited Three-Year Warranty

Your Cube Mate™ Ultrasonic Cool Mist Humidifier is guaranteed to be free of manufacturing defects for a period of three years from the date of purchase under normal use, certain exclusions apply. This limited warranty does not apply to use in any clinical or commercial applications. If the unit fails to operate during the three year period of the original registered purchaser, please call our Customer Care Help Line at 800-622-4714 and we will try to resolve the issue as soon as possible. If Briggs Healthcare® determines that the unit failed to operate due to a manufacturing defect, the unit will be replaced at the option of Briggs Healthcare®. Replacement of the unit is the sole remedy under this limited warranty. This warranty gives you specific legal rights, which vary from state to state.

This Limited Warranty constitutes Briggs Healthcare's only responsibility and obligation to replace the unit. We make no other express or implied warranties, arising by operation of law or otherwise, or any warranty of merchantability or fitness for a particular use or purpose whether or not the use or purpose has been disclosed to Briggs Healthcare® in specifications, drawings or otherwise, and whether or not Briggs Healthcare's products are specifically designed and/or manufactured by Briggs Healthcare® for the buyer's use or purposes, except for the limited warranty stated above. Briggs will not be responsible for any indirect, incidental, special, consequential, or punitive damages or other loss, including, but not limited to, damage to or loss of other property or equipment and personal injuries, whether to purchaser or others. Briggs shall in no event be liable to the purchaser for any amount in excess of the cost of replacement of the unit.

Manufactured for:
Briggs Healthcare
1931 Norman Drive
Waukegan, IL 60085
#91-682-000 06/13
©2013 Briggs Healthcare

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