

Model 15-698-000



**Safety Assured and Comparably Priced
to Glass Thermometers!**

MABIS TinyTemp® Thermometer

One-piece solid plastic construction -
no danger of broken glass

No danger of mercury poisoning

Standard measuring time is 60 seconds

Easy-to-read digital display

Economically priced

Clinical accuracy

Waterproof for easy cleaning

Directions for use

1. Clean before and after use. Unit can be submersed into isopropyl alcohol diluted with water or COLD soapy water or wiped clean. Do not soak for an extended period or sterilize by boiling.
2. Determine measurement site – oral, rectal or under the arm.
3. Press the blue POWER button to activate.
4. The L with the flashing °F in the display window indicates the thermometer is ready for use.
5. When the °F stops flashing, the peak temperature has been reached.
6. Press the blue POWER button to turn the unit off. If left unattended, the unit will automatically shut off after approximately 10 minutes.

IMPORTANT BEFORE USE:

- Normal body temperature is considered 98.6°F (37°C); however variations may occur.
- Measurement times will vary from approximately 1 minute (oral and rectal) up to 10 minutes (under the arm).
- For rectal measurement, lubricate tip and insert a maximum of 1/2" into the rectum.
- Thermometer has a non-replaceable battery with approximately 600 usage hours.
- Thermometer should only be used under the supervision of an adult.
- Do not walk, run or talk during temperature taking.
- Clean the thermometer before and after each use.
- Store the unit in the protective case when not in use.
- Do not bite the probe.
- Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
- Dropping or subjecting your thermometer to strong shocks should be avoided.
- Do not attempt to disassemble the thermometer.
- Thermometer must be safely disposed and replaced when "ERR" displays and measurement is not attainable.

NOTE: Performance of the device may be degraded if: operated outside of stated temperature and humidity range; stored outside of stated temperature and humidity range; thermometer undergoes mechanical shock (drop); patient temperature is below ambient temperature.

Limited One Year Warranty

This MABIS Digital Thermometer is guaranteed to be free of manufacturing defects for a period of one year from the date of purchase under normal use. This warranty applies to the original registered purchaser only. All repairs on your MABIS Digital Thermometer should be sent prepaid to: MABIS Healthcare, Attn: Repair Department, 1931 Norman Drive, Waukegan, IL 60085. MABIS Healthcare shall incur no liability under this warranty for any damage from any cause other than inherent manufacturing defects. ALL IMPLIED WARRANTIES, INCLUDING MERCHANTABILITY AND FITNESS ARE LIMITED TO THE DURATION OF THIS WARRANTY. THIS WARRANTY DOES NOT INCLUDE ANY LIABILITY FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow limitation on how long an implied warranty lasts or do not allow the exclusion or limitation of incidental or consequential damages, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights and you may also have other rights, which vary from state to state.

Please read all instructions before use.

Meets ASTM Standard E 1112

Distributed by:
HealthSmart® International
1931 Norman Drive
Waukegan, IL 60085 USA
Register this product at
LiveHealthSmart.com
©2016 HealthSmart® International

#92-698-000 10/16
Made in China

